


SunCHECK® v4.3 Feature Set

Managing many combinations of Quality Management machines, deliveries and resources can add complexity to day-to-day clinical workflows. SunCHECK Quality Management Platform Version 4.3 helps address these challenges. With over 30 new features and advancements, we are driven to improve your workflows and provide the capabilities you need most – to harness your Quality Management Program.



For more details, see how each new feature supports your team:

Feature	Workflow	Reporting	Analytics	Innovation	Market-Driven
Multi-OIS Support: Multiple OIS systems can connect to a single SunCHECK Instance (On-Premises or Cloud-based)	X			X	
User Permissions Enhancements: Assign users to specific linear accelerators and/or other machines	X				
Transit Dosimetry for Varian Medical Systems® Halcyon™ Treatment Machines: Industry-leading In-vivo dosimetry support for Patient QA				X	X
Secondary Point Dose Calculation: Simple, emergency calculation support for Varian Medical Systems® and Elekta photon Linear Accelerators	X				X
Interrupted Treatment Support: 3D calculation is provided for beam interruption events occurring during delivery	X	X		X	
Asset Management: Manage and link departmental equipment to tasks including TG-51, TRS-398 and TG-142	X				X
External Data Import: Additional data set file types supported (.csv, .xml)	X	X	X		X
Automated Varian Medical Systems® Machine Performance Check (MPC) Data: Import for tracking, trending and reporting within SunCHECK	X	X	X	X	

Varian Medical Systems® is a registered trademark, and Varian™ and Halcyon™ are trademarks of Varian Medical Systems, Inc. Sun Nuclear Corporation is not affiliated with or sponsored by Varian Medical Systems, Inc.



[sunnuclear.com](https://www.sunnuclear.com)

©2024 Mirion Technologies, Inc. or its affiliates. All rights reserved. Sun Nuclear, the Sun Nuclear logo, and other trade names of Sun Nuclear products listed herein are registered trademarks or trademarks of Mirion Technologies, Inc. or its affiliates in the United States and other countries. Third party trademarks mentioned are the property of their respective owners. All Rights Reserved. All data used is best available at time of publication. Data is subject to change without notice.