




## v5.1 Feature Set

With version 5.1 of the SunCHECK® Platform, we are pleased to introduce 15+ new features and advancements, including DoseCHECK™ Monte Carlo and the AdaptCHECK™ module.



See below for an overview of key v5.1 features, and their areas of impact.

Feature	Workflow	Reporting	Analytics	Innovation	User-Driven
<b>Platform Improvements:</b> Continued user interface refresh and tools to validate automation	X				X
<b>Strengthened Security:</b> Updated third-party components and additional penetration testing					X
<b>NEW! Pre-Plan Feasibility (PlanCHECK™):</b> Estimate ability to meet plan objectives prior to creating the treatment plan	X		X	X	
<b>NEW! DoseCHECK™ Monte Carlo:</b> Adds ability to generate 3D Secondary Dose Calculations using a Monte Carlo algorithm for photon and electron energies for linear accelerators, plus Accuray CyberKnife® systems				X	X
<b>Varian Medical Systems® Ethos™ Therapy Integrations:</b>					
<b>PlanCHECK and DoseCHECK Support</b> for initial Plan Checks and Secondary Dose Calculations	X			X	X
<b>PerFRACTION™ Log-based Analysis</b> for Pre-treatment QA and In-vivo Monitoring	X			X	X
<b>NEW! AdaptCHECK™ for Online Adaptive Patient-Specific QA,</b> including Secondary Checks and Log-based Analysis	X	X		X	X
<b>NEW! MLC Log File Analysis</b> during treatment delivery	X	X	X		X

■ Platform-level features
 ■ Patient QA-focused features
 ■ Machine QA-focused features



[sunnuclear.com](http://sunnuclear.com)

©2025 Mirion Technologies, Inc. or its affiliates. All rights reserved. Sun Nuclear, the Sun Nuclear logo, and other trade names of Sun Nuclear products listed herein are registered trademarks or trademarks of Mirion Technologies, Inc. or its affiliates in the United States and other countries. Third party trademarks mentioned are the property of their respective owners. All Rights Reserved. All data used is best available at time of publication. Data is subject to change without notice.

Feature	Workflow	Reporting	Analytics	Innovation	User-Driven
<b>NEW! Excel Data Import</b>	X				X
<b>Surface Guided Radiotherapy Test Support:</b> Record and document SGRT test results for compliance of TG-302, TG-147 and ESTRO-ACROP standards	X	X		X	
<b>Additional Phantom Support:</b> Mobius MC2, Quart DVT_VN included parameters		X			X
<b>Image-based QA Task Support for Varian Medical Systems® Halcyon™ and Ethos™ Systems:</b> Acquire, analyze and trend leaf speed and dose rate/gantry rotation	X			X	X

■ Platform-level features    
 ■ Patient QA-focused features    
 ■ Machine QA-focused features



Learn more about  
 SunCHECK v5.1 on the  
**Release Highlights**  
 page at [sunnuclear.com](http://sunnuclear.com)

*Varian Medical Systems® is a registered trademark, and Varian™ and Halcyon™ are trademarks of Varian Medical Systems, Inc. Sun Nuclear Corporation is not affiliated with or sponsored by Varian Medical Systems, Inc.*