




v4.2.1 Feature Set

This release of SunCHECK is a release for both On-Premise and SaaS/Cloud users of the SunCHECK Platform. SunCHECK v4.2.1 delivers a feature-rich release with multiple key new capabilities and 100+ enhancements, to improve workflows. These new features and user-requested improvements set the stage for clinics to optimize their Quality Management programs.



See below for an overview of key v4.2.1 features, and their areas of impact.

Feature	Workflow	Reporting	Analytics	Innovation	User-Driven
Over 100 valuable user-requested enhancements and reported issues addressed	X	X			X
Single Sign-On (SSO) Support: Manage authentication for SunCHECK access through your facility's Microsoft Active Directory solutions	X			X	
Data Management Improvements: Updated architecture to rely on file system storage	X				
EPID-based 2D pre-treatment and relative in-vivo dosimetry for Varian Medical Systems® Halcyon™ System: Industry leading In-vivo dosimetry support for Patient QA				X	X
Multiple CT-to-Electron Density Table support: Improved selection of multiple CT or CBCT protocols for dose calculations	X				X
Template Sharing: Improvements to allow sharing of templates with other machines	X	X			X

■ Platform-level features
 ■ Patient QA-focused features
 ■ Machine QA-focused features

Please note: SunCHECK Version 4.2.1 is required to be able to update to Version 5.0.

Varian Medical Systems® is a registered trademark, and Varian™ and Halcyon™ are trademarks of Varian Medical Systems, Inc. Sun Nuclear Corporation is not affiliated with or sponsored by Varian Medical Systems, Inc.



sunnuclear.com

©2024 Mirion Technologies, Inc. or its affiliates. All rights reserved. Sun Nuclear, the Sun Nuclear logo, and other trade names of Sun Nuclear products listed herein are registered trademarks or trademarks of Mirion Technologies, Inc. or its affiliates in the United States and other countries. Third party trademarks mentioned are the property of their respective owners. All Rights Reserved. All data used is best available at time of publication. Data is subject to change without notice.